



Environment Friendly



Energy Saving (VSD)



Quiet



DENTAL SUCTION SYSTEM

- Flexible installation, easy maintenance.
- Low running noise, suitable for dental clinic.
- Stainless made boosting pump can provide longer product life.
- Recycle water cooling system saving energy and environment friendly.
- Strong and sensitive suction, providing sufficient negative pressure suction at any time.
- Equipped with a delayed shutdown design, reducing the frequency of motor restarts to extend motor life.
- Environmental protection without consumables, regular cleaning of filters and pipes can maintain suction.



PM7204



PM7204S

Model	Cooling	Voltage	kw	MAX VACUUM mmHg / Kpa	dBA	kg	cm	Chair / Dentist	Remark
PM7204	Water	220V	1.10	150 / \leq -20	55	27	36x30x40	5 / 4	
PM7204S	Water	220V	1.10	150 / \leq -20	52	46	45x36x53	5 / 4	

PUMA INDUSTRIAL CO., LTD.

OFFICE / FACTORY :

No. 427, Sec. 2, Doujhong Rd., Tianjhong Township, Changhua County 520, Taiwan

TEL:886-4-8761811

FAX:886-4-8752211

e-mail:sales@pumabiomedical.com

www.pumabiomedical.com



ISO-9001 / 13485 CE FDA REACH

(SPECIFICATION SUBJECT TO CHANGE WITHOUT NOTICE.)

PMD-1930L-11D
© 2019