



VACUUM PUMP

- Instant react speed.
- Easy maintenance.
- Flexible installation.
- Powerful.



VS01



VS08

									
Model	Voltage	kw	ℓ / min	mmHg / Kpa	dBA	kg	cm	Chair / Dentist	Remark
VS01	110 / 220V	0.07	22	680 / ≤-90	58	4	28x16x25	2 / 1	
VS08	110 / 220V	0.56	107	700 / ≤-93	61	10	30x22x23	3 / 2	

PUMA INDUSTRIAL CO., LTD.

● OFFICE / FACTORY :

No. 427, Sec. 2, Doujiong Rd., Tianjiong Township, Changhua County 520, Taiwan

TEL:886-4-8761811

FAX:886-4-8752211

e-mail:sales@pumabiomedical.com

www.pumabiomedical.com

