



牙科專業無油空壓機

- 高效率，低噪音。
- 自動運轉，節約能源。
- 專利設計，堅固耐用。
- 適用於牙科診所。



MG130V



LG230V



WD230V



PF330V



WD460 / 8110



OF250VA

型號	電壓 / 相數 Volt / Ph	馬力 / HP	氣缸數 NO.	排氣量 l/min @ 3kg/cm ² 50Hz / 60Hz	使用壓力 Kg/cm ² / psi	儲氣桶容量 l / Gal.	噪音值 dBA	淨重 kg / Lbs.	尺寸 cm	備註 請洽業務 / 牙科使用 Chair / Dentist
MG130V	220V / 1	0.38 / 0.5	1	46 / 57	8 / 115	30 / 8	73	31 / 68	46x46x66	2 / 1
LG230V	220V / 1	1.10 / 1.5	1	82 / 90	8 / 115	30 / 8	74	38 / 83	46x46x76	3 / 2
PF330V	220V / 1	1.50 / 2.0	2	149 / 172	8 / 115	30 / 8	78	53 / 117	46x46x76	4 / 3
OF250VA	220V / 1	1.84 / 2.5	2	128 / 153	8 / 115	50 / 13	81	46 / 101	50x46x85	4 / 3
WD230V	220V / 1	1.10 / 1.5	2	101 / 122	8 / 115	30 / 8	63	37 / 82	42x39x63	3 / 2
WD460	220V / 1	2.20 / 3.0	4	187 / 229	8 / 115	58 / 15	65	65 / 143	96x37x65	6 / 4
WD8110	220V / 1	4.40 / 6.0	8	374 / 458	8 / 115	110 / 30	66	118 / 260	123x43x73	10 / 8

● 另有無油式機體可供選購，有任何需求請與我們聯絡。

合正機械股份有限公司醫器廠 **PUMA INDUSTRIAL CO., LTD.**

● 公司 / 工廠

52046 彰化縣田中鎮斗中路二段427號3樓

TEL: 886-4-8761811

FAX: 886-4-8752211

e-mail: sales@pumabiomedical.com

www.pumabiomedical.com

